

BRONCOFLEX® AGILE

Single-use sterile video-bronchoscope



axessvision
part of The Surgical Company

HIGH QUALITY VISION & SAFETY.
ANYTIME, ANYWHERE.

BRNCOFLEX® AGILE

GUIDING AND INSPECTION

Used in combination with endotracheal tubes, Broncoflex® Agile makes it easier to guide tubes through the anatomical structures and insert them thanks to the rigidity of its braided Pebax® insertion tube.

The correct positioning of double-lumen endotracheal tubes can be verified for sizes of 35 Fr or more.



ERGONOMICS & MANOEUVRABILITY

The precise, flexible deflection mechanism offers up and down distal angulation of 220°, facilitating bronchial exploration, particularly of the right upper lobe of the lung.



HIGH-DEFINITION SYSTEM

The whole system, from the camera to the visualisation platform, contributes to the display of a high-definition image with, in particular, an algorithm especially developed for detailed contours. Brightness management enables good near and distal vision.



SECRETION MANAGEMENT

The unique profile of the distal tip, with its wide opening, helps facilitate the suction of secretions and fluids.

TECHNICAL CHARACTERISTICS

| | | |
|--------------------------|---|----------------------------------|
| OPTICS | Field of vision direction | 0° |
| | Field of vision | 87.5° |
| | Field depth | 5-50mm |
| | Image resolution (pixels) | 400x400 |
| | Lighting system | 2 LEDs |
| INSERTION TUBE | High/low deflection angle | 220°/220° |
| | External diameter of insertion tube | 3,6mm |
| | Distal end outer diameter | 3,9mm |
| | Working length | 605mm |
| WORKING CHANNEL | Internal diameter | 1,4mm |
| STERILISATION | Ethylene oxide | |
| INSTRUMENT COMPATIBILITY | Compatible with intubation tube (internal diameter) | 5mm |
| | Double-lumen tube (internal diameter) | 35Fr |
| REFERENCES | Broncoflex Agile Broncoflex Vortex Screeni HD | 20030001 10030001 30030000 |